

## Product Information

AmnioPrime

Complete Medium for Cultivation of Amnion and Chorionic Villi Cells

Cat. No.: APR-B (100 ml)



### General Information

AmnioPrime is a ready to use medium for in vitro diagnostic use. It has been specifically developed for the cultivation of human primary amnion and chorionic villi cells, which are intended for the preparation of karyograms, fluorescence *in situ* hybridization and other cytogenetic methods.

The medium is supplied frozen and contains preselected serum, L-glutamine and antibiotics.

Appearance	Clear yellow to red frozen liquid
CO <sub>2</sub> concentration, optimum	5%
Storage and shelf life	Store at ≤-15°C protected from light. Once opened, store at +2°C to +8°C and use within 2 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	+37°C in water bath and swirl gently to homogenize.

### Formulation

The Formulation is based on the basal medium MEM Alpha Modification. Note that the liquid proprietary formulation of the medium already contains preselected fetal bovine serum, L-glutamine and 50 µg/ml gentamycin.

### Instructions for Use

The medium may be used in both open and closed culture systems.

#### *Important information*

Supplementation of AmnioPrime Medium is neither necessary nor recommended. It is recommended to use cells from 2.5 ml of amniotic fluid per one coverslip. The following protocol and the volumes indicated are only general guidelines for use. This high quality medium can be used within established procedures. It is up to the user to adopt either parts or all of the optimized protocol described below.

#### *In situ Culture of Amniotic Fluid Cells*

1. Concentrate the cells by centrifugation of the amniotic fluid: Centrifuge 20 ml of amniotic fluid at 750 rpm for 10 minutes.
2. Carefully decant the amniotic fluid from the cell pellet into a sterile test tube.
3. Resuspend the cell pellet with 2 ml of amniotic fluid.
4. Add 2 ml of AmnioPrime Medium and swirl gently.
5. Culture 0.5 ml of the cell suspension on each coverslip in a tissue culture dish.
6. Incubate cultures at +37°C in a 5% CO<sub>2</sub> atmosphere.
7. Add 2 ml of AmnioPrime Medium to each culture on day 2.
8. Check cultures for growth after 4 to 5 days. Feed cultures once growth has been observed. To feed cultures, carefully aspirate all of the exhausted culture medium and replace with 2 ml of fresh AmnioPrime Medium. Recommendation: feed cultures every 2 days.
9. When the cultures have colonies of sufficient size, proceed with harvesting.

## Product Information

AmnioPrime

Complete Medium for Cultivation of Amnion and Chorionic Villi Cells

Cat. No.: APR-B (100 ml)



10. For best results, feed cultures with AmnioPrime Medium the day before the harvest.

### Flask Method Culture of Amniotic Fluid Cells

1. Use the same procedure as for the in situ culture, with the following adaptations:
2. Re-suspend the cell pellet with 4 ml of amniotic fluid. Add 16 ml of AmnioPrime Medium and swirl gently.
3. Culture 5 ml per T25 flask. Place the cap loosely on the flask and incubate undisturbed at +37°C in a 5% CO<sub>2</sub> atmosphere.
4. Check all flasks for growth after 5 days.
5. For best results, feed cultures with AmnioPrime Medium the day before the harvest.

### Recommendations for Closed Systems

AmnioPrime Medium may be used in closed culture systems as long as the physiological pH of 6.9 to 7.4 is maintained. Closed systems depend on adequate buffering capacity of media.

Method 1: Supplement AmnioPrime Medium with 2% (v/v) sterile 1.0 HEPES solution. The HEPES solution must be set to pH 7.0 at +20°C. HEPES supplemented medium can subsequently be used on cells in closed culturing flasks.

Method 2: Pre-equilibrate the flask containing AmnioPrime Medium and cells at +37°C in a 5% CO<sub>2</sub> atmosphere for 1 hour prior to closing the flask.

Method 3: Flush each culture flask containing AmnioPrime Medium and cells with 5% CO<sub>2</sub> – 95% air through 0.2 µm sterile filter for 20 seconds. Tighten the caps and incubate the flasks at +37°C.

### Precautions and Disclaimer

For *in vitro* diagnostic use.

The medium is not intended for therapeutic use.

Do not use if a visible precipitate is observed in the medium.

Use of AmnioPrime Medium does not guarantee the successful outcome of any prenatal diagnostic testing.

Do not use AmnioPrime Medium beyond the expiration date indicated on the product label.

### Help needed?

If you have any further questions regarding this product please do not hesitate to contact our cell culture experts by email ([techservice@capricorn-scientific.com](mailto:techservice@capricorn-scientific.com)) or phone (+49 6424 944640).