

## Product Information

FreezeMe ZERO

DMSO-free Cryopreservation Medium for Cell Lines

sterile-filtered

Cat. No. FMZ-50ML (50 ml), 5-FMZ-10ML (5 x 10 ml), FMZ-10ML (10 ml)

### General Information

FreezeMe ZERO is a fully chemically defined, DMSO-free cryopreservation medium developed for reliable long-term storage of mammalian cell lines. The formulation is specifically designed to preserve cellular integrity during freeze-thaw cycles while avoiding transcriptomic, metabolic, and epigenetic perturbations associated with dimethyl sulfoxide (DMSO).

Unlike conventional cryopreservation media, FreezeMe ZERO eliminates DMSO entirely, thereby minimizing freeze-thaw-induced molecular alterations that may impact long-term cell line stability and performance. The medium supports high-density cryopreservation and enables robust post-thaw recovery, making it particularly suitable for the generation and storage of Master Cell Banks (MCBs).

### Applications

- Cryopreservation of mammalian cell lines for Master and Working Cell Banks
- Long-term storage of production-relevant cell lines
- High-density cell banking workflows
- DMSO-free preservation strategies for sensitive cell lines

### Product Specifications

Appearance	Clear solution
Specifications	<ul style="list-style-type: none"> <li>- DMSO-free</li> <li>- Chemically defined</li> <li>- Serum-free</li> <li>- Animal derived component-free</li> </ul>
Storage and Shelf Life	Store at +2°C to +8°C Protect from light.
Shipping Conditions	Ambient

### Instructions for Use

#### General Notes:

- FreezeMe ZERO is suitable for adherent and suspension cell lines. For freezing of adherent cell lines, detach the cells with a gentle dissolving reagent, such as Accutase (ACC-1B) and resuspend the cells in fresh culture medium
- Cell viability should be ideally > 90% before cryopreservation

#### Cryopreservation of cells:

1. Collect the cells in a sterile tube and centrifuge for 5 min at 400 x g
2. Carefully remove the supernatant, resuspend in DPBS buffer (e.g. PBS-1A) and determine the cell count
3. Pellet for 5 min at 400 x g and resuspend in FreezeMe ZERO to reach a cell density of  $1 \times 10^6$  to  $1 \times 10^8$  cells/ml.

## Product Information

FreezeMe ZERO

DMSO-free Cryopreservation Medium for Cell Lines

sterile-filtered

Cat. No. FMZ-50ML (50 ml), 5-FMZ-10ML (5 x 10 ml), FMZ-10ML (10 ml)

4. Aliquot the cell suspension to a freezing volume of 0.5 – 1.5 ml. If freezing in larger volumes such as bags is intended, please use programmed cooling curves
5. Place the aliquots to a -80°C freezer. After 12 hours, transfer the aliquots to liquid nitrogen for long-term storage

### Thawing of cryopreserved cells:

1. Immediately, after removing from liquid nitrogen, place the cryopreserved cells in a 37°C water bath and thaw cells rapidly for 2 min
2. Dilute the cells (1:10) to a conical tube with fresh, prewarmed culture medium
3. Invert 10 times to equilibrate
4. Centrifuge for 5 min at 400 x g and remove supernatant
5. Gently resuspend the cells in fresh medium and determine viability

### Precautions and Disclaimer

---

This product is for research and further manufacturing use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### Help Needed?

---

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).