

## Product Information

G-418 Sulfate Solution, 100 % Activity (50 mg/ml)

Cat. No. G418-H (10 ml), G418-B (100 ml), G418-10X (10x10 ml)

### General Information

G-418 is used in the selection and maintenance of eukaryotic cells, stable transfected with neomycin resistance genes. G-418 is an aminoglycoside antibiotic, related to Gentamicin, and exhibits toxicity towards both eukaryotic and prokaryotic cells. It is produced by *Micromonospora rhodorangea* and acts by binding the ribosome, thus inhibiting protein synthesis in both prokaryotic and eukaryotic cells.

### Product Specifications

Appearance	Clear frozen liquid
CAS No.	108321-42-2
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at $+4^{\circ}\text{C}$ and use within 4-6 weeks.
Shipping conditions	Frozen (Dry Ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$ . Swirl gently to homogenize.
Working concentration	Recommended final concentration (0.1 – 1.0 mg/ml) depending on the cell type: <ul style="list-style-type: none"> <li>• HeLa: 200 – 600 <math>\mu\text{g}/\text{ml}</math></li> <li>• 3T3 cells: 500 – 1000 <math>\mu\text{g}/\text{ml}</math></li> <li>• CHO: 200 – 400 <math>\mu\text{g}/\text{ml}</math></li> <li>• HEK 293: 500 – 800 <math>\mu\text{g}/\text{ml}</math></li> <li>• Jurkat cells: 600 – 700 <math>\mu\text{g}/\text{ml}</math></li> </ul>

### Instructions for Use

- Do not use G-418 with antibiotic/antifungal preparations (e.g. Pen/Strep). These agents are competitive inhibitors of G-418. Other antibiotics are potentially cross-reactive as well.
- Good laboratory practice requires optimal concentration of biologically active G-418 to select and maintain cells. This must be determined for each set of growth conditions. G-418 is used in the concentration range of 100 – 200  $\mu\text{g}/\text{ml}$  for bacteria, or 200 – 500  $\mu\text{g}/\text{ml}$  for most mammalian cells. Concentrations of G-418 required for maintenance of selected cell lines are typically  $\leq 50\%$  compared to selection.
- It is recommended, whenever experimental conditions are altered, the optimal concentration of the product should be re-evaluated.

### Precautions and Disclaimer

This product is for research use and further manufacturing only.

### Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).