

Product Information

Water, suitable for cell culture, sterile-filtered
Cat. No. WAT-1000ML (1000 ml), WAT-500ML (500 ml)

General Information

High quality water for use in cell culture applications, purified by distillation, ion exchange via reverse osmosis, deionization, and ultraviolet disinfection. It can be used for a variety of applications, such as preparation of cell culture media, buffers, and laboratory reagents, rinsing bioreactors, plasticware and glassware, or for methods like chromatography (HPLC) and spectrometry.

Product Specifications

Appearance	Clear colorless liquid
Storage and shelf life	Store at +2°C to +25°C upon arrival. Once opened store at +2°C to +25°C and use within 24 months.
Shipping conditions	Ambient

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).