

Product Information

CHO|ONE A, Adaptation Medium, w/o Insulin, w/o L-Glutamine, with Pluronic™
Cat. No. CHOA-500ML (500 ml)

General Information

CHO|ONE A allows the cultivation or adaptation of recombinant CHO cells in small-scale shaker bottles or spinner flasks. Lacking certain metabolites, CHO|ONE A exerts a strong selection pressure on the cells. In addition, its formulation simplifies and shortens the time-consuming adaptation to fed-batch culture.

Product Specifications

Appearance	Clear yellow orange solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).