

## Product Information

CHO|ONE Feed 1B

Feeding Supplement for CHO Cells with Pluronic™, w/o Insulin, w/o L-Glutamine

Cat. No. CHOF1B-250ML (250 ml)

### General Information

---

CHO|ONE Feed 1B is a feeding solution with modified salt composition that can be universally combined with all CHO|ONE Media formulations in a fed batch process.

This feed formulation provides a selection of nutrients with optimized salt concentration to facilitate the control of osmolality and metabolite accumulation during fed-batch culture. It provides high cell viability and productivity, especially for CHO clones sensitive to osmotic alterations. All CHO|ONE Feed 1 variants can be combined with CHO|ONE Feed 2 for an efficient feeding strategy in fed-batch processes.

### Product Specifications

---

Appearance	Clear red to orange solution
Storage and Shelf Life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping Conditions	Ambient

### Formulation

---

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

### Precautions and Disclaimer

---

This product is for research use and further manufacturing only.

### Help Needed?

---

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([techservice@capricorn-scientific.com](mailto:techservice@capricorn-scientific.com)) or phone (+49 6424 944640).