

Product Information

CHO|ONE Feed 1C

Feeding Supplement for CHO Cells with Pluronic™, w/o Insulin, w/o L-Glutamine

Cat. No. CHOF1C-250ML (250 ml)

General Information

CHO|ONE Feed 1C is a low-nutrient feeding solution that can be universally combined with all CHO|ONE Media formulations in a fed batch process.

This light, metabolically gentle feed formulation provides a reduced nutrient load. Therefore CHO|ONE Feed 1C helps to minimize metabolic stress or overload, for sensitive CHO clones. It can be used to identify minimal feeding requirements while maintaining cell viability. All CHO|ONE Feed 1 variants can be combined with CHO|ONE Feed 2 for an efficient feeding strategy in fed-batch processes.

Product Specifications

Appearance	Clear red to orange solution
Storage and Shelf Life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping Conditions	Ambient

Formulation

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).