

## Product Information

### CHO|ONE L

Expression Medium with Pluronic™, w/o HT, w/o Insulin, w/o L-Glutamine

Cat. No. CHOL-500ML (500 ml)

### General Information

---

CHO|ONE L allows the cultivation of recombinant CHO cells for fed-batch cultures in spinner flasks or large-scale bioreactors. This light nutrient formulation with reduced amino acid and glucose levels is ideal for identifying minimal requirements, reducing metabolic stress, or screening nutrient sensitivities.

CHO|ONE L already leads to increased cell growth and protein yield, but combined with CHO|ONE Feeding supplements the full potential of the Media System can be deployed.

### Product Specifications

---

Appearance	Clear yellow orange solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

### Formulation

---

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

### Precautions and Disclaimer

---

This product is for research use and further manufacturing only.

### Help Needed?

---

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([techservice@capricorn-scientific.com](mailto:techservice@capricorn-scientific.com)) or phone (+49 6424 944640).