

Product Information

CHO|ONE Feed 2
Supplement for CHO cells, w/o Insulin, w/o L-Glutamine
Cat. No. CHOF2-10ML (10 ml), CHOF2-50ML (50 ml)

General Information

Capricorn's CHO|ONE Feed 1 and Feed 2 are feed additives especially matched to the core medium CHO|ONE E. Combining those increases fed batch productivity 3- to 5-fold compared to batch culture conditions. Proposed protocols facilitate easy and efficient identification of a feeding strategy. CHO|ONE Feed 1 contains concentrated amino acids, vitamins, salts, trace elements and a carbon source. CHO|ONE Feed 2 is a balanced formulation of selected amino acids in a concentrated form.

Product Specifications

Appearance	Clear solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).