

Product Information

Dulbecco's Modified Eagle Medium (DMEM) High Glucose (4.5 g/L), with L-Glutamine, with Sodium Pyruvate
 Cat. No. DMEM-HPA (500 ml)

Product Specifications

Appearance	Clear red orange solution
CO ₂ concentration, optimum	8.5 %
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/L
Amino Acids:	
Glycine	30.00
L-Arginine HCl	84.00
L-Cystine 2 HCl	62.57
L-Glutamine	584.00
L-Histidine HCl H ₂ O	42.00
L-Isoleucine	105.00
L-Leucine	105.00
L-Lysine HCl	146.00
L-Methionine	30.00
L-Phenylalanine	66.00
L-Serine	42.00
L-Threonine	95.00
L-Tryptophan	16.00
L-Tyrosine 2 Na 2H ₂ O	103.79
L-Valine	94.00

Components	Concentration mg/L
Vitamins:	
Choline chloride	4.00
D-Calcium Pantothenate	4.00
Folic Acid	4.00
myo-Inositol	7.20
Nicotinamide	4.00
Pyridoxal HCl	4.00
Riboflavin	0.40
Thiamine HCl	4.00
Inorganic Salts:	
CaCl ₂ 2 H ₂ O	265.00
Fe(NO ₃) ₃ 9 H ₂ O	0.10
KCl	400.00
MgSO ₄ 7H ₂ O	200.10
NaCl	6400.00
NaHCO ₃	3700.00
NaH ₂ PO ₄ 2H ₂ O	141.73
Other Components:	
D-Glucose	4500.00
Phenol Red Sodium Salt	15.90
Sodium Pyruvate	110.00

Precautions and Disclaimer

This product is for research use only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).