

Product Information

Other Animal Sera

Cat. No. GOA-1B, GOA-1A, POR-1B, POR-1A, MOU-1B, RAT-1B, RAB-1B

General Information

Other animal sera are used in various applications, e.g. autologous cell culture systems or immunohistological methods.

Porcine serum is often used in cultures that proliferate mycoplasma or viruses. Goat and Sheep Serum are used in immunohistochemistry and to eliminate non-specific binding reactions.

Applications:

- Autologous cell culture systems
- Immunohistochemistry
- Immunoblotting
- Immunofluorescence
- Immunoprecipitation

Test Parameter:

- pH
- Osmolality
- Endotoxin
- Total protein
- Hemoglobin
- Mycoplasma
- Sterility

Product	Volume	Cat. No.	Storage
Goat Serum	100 ml 500 ml	GOA-1B GOA-1A	Store at $\leq -15^{\circ}\text{C}$.
Porcine Serum	100 ml 500 ml	POR-1B POR-1A	Store at $\leq -15^{\circ}\text{C}$.
Mouse Serum	100 ml	MOU-1B	Store at $\leq -15^{\circ}\text{C}$.
Rat Serum	100 ml	RAT-1B	Store at $\leq -15^{\circ}\text{C}$.
Rabbit Serum	100 ml	RAB-1B	Store at $\leq -15^{\circ}\text{C}$.

Precautions and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).