

Product Information

Fetal Bovine Serum Low Endotoxin, Collected in South America
Cat. No. FBS-LE-12A (500 ml), FBS-LE-12B (100 ml)

General Information

Serum is the blood component obtained after coagulation and removal of cellular components. Besides serum proteins it contains growth factors, amino acids and hormones, etc. This comprehensive mix makes serum one of the most important supplements for supporting cell growth and proliferation in *in vitro* cell culture. Of special interest is fetal bovine serum (FBS). It is especially rich in growth factors and low in antibodies, which may influence the cell culture work.

Fetal Bovine Serum Low Endotoxin is a selected FBS featuring an endotoxin level lower than 1 EU/ml, making it especially suitable for sensitive cells, primary cell culture, and hybridoma cloning.

Product Specifications

Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at $+4^{\circ}\text{C}$ and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Quality Control

Only serum batches which pass our strict quality control are released for sale. Standard parameters which are determined include pH, osmolality, content of protein, albumin, IgG and hemoglobin, endotoxin level, sterility, mycoplasma detection and virus testing.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).