

Product Information

Rabbit Serum
sterile-filtered
Cat. No. RAB-1A (500 ml), RAB-1B (100ml)

General Information

Rabbit Serum is a versatile supplement for various applications in cell culture, immunoassays, and microbiological research. Derived from healthy New Zealand rabbits, it provides essential growth factors, proteins, lipids, and hormones that support cell attachment, proliferation, and overall viability. Rabbit Serum can serve as an alternative to fetal bovine serum (FBS) for selected cell lines and specialized assays.

Product Specifications

Appearance	Amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at $+4^{\circ}\text{C}$ and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Precautions and Disclaimer

This product is for research use and further manufacturing only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).